

## Research Order Form

- -Complete the form below for <u>each separate</u> delivery date and email to <u>Inquiry@BioSharingNetwork.com</u>.
- -Please submit orders at least 5 business days prior to the time products will be shipped or picked up.
- -Ship/pick-up times are **4:00pm** unless otherwise stated in the form below.
- -Orders submitted 3 or 4 days before ship/pickup time will incur a \$25 stat fee.
- -Orders submitted 2 days before ship/pickup time will incur a \$40 stat fee.
- -Orders received 1 day before ship/pickup date **cannot** be filled.
- -For buffy coat and platelet orders on days with low predicted whole blood draws, you may be asked to choose a different ship date or different product.

Re	equesting Company:											
Primary (	Contact and/or Lab Name:											
Pe	erson Placing Order:											
	Delivery Date:											
	Product Use:	Research Further Manufa	acturing	Other								
Catalog Number		Description	Product Price	Donor Demo- graphics	Quantity							
BUFFY COATS - Whole blood-derived, manufactured <8 hours of collection; 1-2.5x10° leukocytes (30-80% PMNs) in ~50 mL plasma												
BC-8L	Buffy Coat; available at 4:00	\$100	<b>\$10</b>									
BC-8E	Buffy Coat; available at 10:00	Dam day <u>after</u> collection	\$100	<b>\$10</b>								
LEUKOREDUCTION SYSTEM (LRS) CHAMBERS - Byproduct of platelet apheresis; concentrated PBMCs (1-2.5x10° PBMCs, <10% PMNs) in ~10 mL blood; "buffy cone"  ** Available only Tuesday – Friday after 12:00pm. Cell counts cannot be guaranteed when specific donor demographics are requested.												
LRSC-01	LRS chamber, sterilely sealed	; available at 4:00pm day of collection	\$80	<u></u> \$10								
LRSC-02	LRS chamber, sterilely sealed	\$65	Included									
LRSC-03	LRS chamber; drained, cells s	\$90	Included									
ISOLATED PBMCs - Isolated by Ficoll from LRS chambers the day after collection  ** Available only Tuesday – Friday after 12:00pm. Cell counts cannot be guaranteed when specific donor demographics are requested.												
PBMC-01	PBMCs suspended in buffer,	\$200	Included									
PBMC-02	PBMCs, cryopreserved; 108/v	\$125	Included									
PLATELETS - Whole blood-derived platelet rich plasma (PRP); CPD anticoagulant												
PRP-8S	Single-donor "random" PRP;	avg of ~6x10 <sup>10</sup> platelets in 40-60 mL plasma	\$100	<u></u> \$10								
PRP-8C	Platelet count provided on si	ngle-donor PRP unit	\$25	N/A								
PRP-24P	Pooled PRP from 5-6 donors;	avg of ~3x10 <sup>11</sup> platelets in 250-350 mL plasma	\$200	N/A								
PLASMA - Whole bloc	od-derived, 250-300 mL/unit; C	PD anticoagulant										
RPFF	Plasma isolated from whole b	\$60	<b>\$10</b>									
RED BLOOD CELL (RBC) AND WHOLE BLOOD UNITS - (subject to availability)												
RBC-01	Packed RBC unit, fresh, Adso	l anticoagulant; available day after collection	\$280	<u></u> \$10								
RBC-02	Packed RBC unit, expired; Ac	dsol anticoagulant	\$25	<u></u> \$10								
WB-01	Whole blood unit, fresh; CPD	anticoagulant; available day of collection	\$550	<b>\$10</b>								

1

RESIDUAL SAMPLE TUBES - Plasma (EDTA, purple top) or Serum (Red top), 48-72 hours after collection, not sterile											
EDTA-01	, , , ,	DTA anticoagulant	\$20	\$10							
EDTA-02	Two 6 mL EDTA purple-top tubes from same donor						\$10				
SRT-01	Single 6 mL red-to	p tube containing clot	\$30	<u></u> \$10							
SRT-02	Two 6 mL red top	wo 6 mL red top tubes from the <u>same donor</u>									
SRT-03							Included				
POOLED TUBE SERUM, MULTI-DONOR - 450-500 donors/L of serum, pooled 72-96 hours from collection, verified negative infectious disease test results											
PTS-01	3.5 L Pooled tube serum, not sterile, frozen					\$1592	N/A				
PTS-02	1.0 L Pooled tube	serum, not sterile, froze	\$455	N/A							
PTS-03	1.0 L Pooled tube	serum, not sterile, refrig	\$455	N/A							
PTS-04	1.0 L Pooled tube	0 L Pooled tube serum, filter sterilized, refrigerated (made on-demand)									
PTS-05	1.0 L Pooled tube serum, filter sterilized, frozen						N/A				
PTS-06	500 mL Pooled tube serum, filter sterilized, frozen						N/A				
PTS-07	250 mL Pooled tube serum, filter sterilized, frozen						N/A				
PTS-08	50 mL Pooled tube serum, filter sterilized, frozen						N/A				
DONOR SEL	ECTION AND RECR	UITMENT –									
Sel-01	Age:G	n already-collected produc ender: Race	ə:	CMV	:	Quote	Included				
Rec-01	Rec-01 For the recruitment of a donor for a product collection based on specific demographics  Age: Gender: Race: CMV:  Other:  COMMENTS OR SPECIAL INSTRUCTIONS:						Included				
COMMENTS OR SPECIAL INSTRUCTIONS:											
		Shippin	g/Pick-	up Info	rmation						
Product(s) will be picked up at BSN Center. Phone number to call/text when order is ready for pick up:  For pick-up orders, please fill out all fields below except FedEx Account Number and Delivery Time.											
Note: A ship	ping charge will be a	ıdded if no FedEx accoun	t number	is provide	ed						
FedEx Acc	count Number:										
Shipping Temperature:		Ambient	Ambient		Refrigerated		☐ Frozen				
Delivery Time:						10:30am					
Ship-to Company:											
Shipping A	. ,										
City: State:		State:	Zip:								
Contact Person:		01010.		Phone Number:							
Billing Information - PO or Credit Card (Visa/MasterCard only)											
PO Number:											
CC Name:				ımber:							
Exp. Date:			CVVC								